

FALL PREVENTION PROGRAM

VIRTUAL BED RAILS

Fall Prevention and Patient Safety & Security is our Priority

NurseView, our continual monitoring system, provides a level of visual contact between nurse and patient. In the event of a potential incident, nursing staff can anticipate a patient action or need and get to the bedside in seconds as opposed to minutes or even hours.


If a patient is identified as high-risk for falls, nursing staff will also activate Virtual Bed Rails as added protection. **Virtual Bed Rails** are noninvasive motion sensitive borders that trigger an onscreen audible alarm when a

defined movement crosses the bed rail to warn the nursing staff of a potential fall that could result in an injury. As an added security feature, room images can be recorded and archived for review and to support legal mitigation.

CareView is the next generation of patient care monitoring, using innovative technology and solutions to address the prevention of falls, patient safety & security, and to provide a more rewarding and pleasant hospital experience.

Clinical Studies using CareView's Virtual Bed Rail Product have resulted in fall reduction in excess of 33%.

Virtual Bed Rails provides *smart* zone coverage specific to a patient's normal physical activity while in hospital care.



CareView's touch-screen monitor provides visual contact from the Nurse's Station and reduces unnecessary steps.



0.176	0	0.278	0	0.193
0	0	0.046	0	0

CAREVIEW
COMMUNICATIONS

CONNECTING PATIENTS, FAMILIES, AND
HEALTHCARE PROVIDERS


BEDSORE MANAGEMENT PROGRAM

ULCER PREVENTION PRODUCT

The development of pressure ulcers results in additional health concerns, longer stays, and increased medical cost at the expense of the hospital.

The CareView **Bedsore Management Program** ensures that patients who are at risk for developing pressure ulcers are repositioned, in the event that the patient is unable to move on their own in a timely manner. If a patient has not been moved, an on-screen visual and audible alarm is triggered at the nurse's station monitor prompting the caregiver to change the patient's position on a regular recurring prescribed schedule.

Providing a more rewarding and pleasant hospital experience using innovative technology and solutions to address the prevention of bedsores and patient safety & security is what makes CareView the **next generation of patient care monitoring.**



A histogram prompts an onscreen audible alarm to indicate the patient has not moved in a prescribed amount of time.

CARE VIEW
COMMUNICATIONS

Phone: 972.943.6050

WWW.CARE-VIEW.COM

405 State Highway 121 Bypass, Suite B-240 • Lewisville, TX 75067 • Fax: 972.403.7659 • info@Care-View.com